

User Manual for KetoSmart Device

Introduction:

Welcome to the user manual for our Acetone Monitoring Device. This manual is designed to guide you through the proper usage and maintenance of the device to ensure accurate and safe monitoring of acetone levels in your breath.

Device Description:

Our Acetone Monitoring Device is a portable device designed to measure acetone levels in your breath. It consists of a handheld unit with three indicator lights: green for normal acetone levels, orange for moderate acetone levels, and red for high acetone levels. The device also contains GPS functionality to detect your location during use.

Getting Started:

To begin using the device, ensure it is fully charged using the provided charging cable. Once charged, power on the device by pressing the power button.

Operation:

Exhale into the device: blow into the mouthpiece.

Wait for the indicator lights: the LED will illuminate for 1 minute during the test.

- Green light indicates normal acetone levels.
- Orange light indicates moderate acetone levels.
- Red light indicates high acetone levels.

Check your mobile app: you will receive a notification on your mobile app for more detailed information.

Repeat the test: after the LED turns off, you can blow into the device again.

Location Detection: for location detection during the test, remove the cap from the device before exhaling to allow the antenna to detect your location.

Maintenance and Care:

Clean the mouthpiece of the device after each use using a mild detergent and water.

Store the device in a cool, dry place when not in use to prevent damage.

Safety Information:

Only use the device as instructed in this manual.

Keep the device away from moisture and extreme temperatures.

Do not attempt to disassemble or repair the device yourself.

Consultation with a Doctor:

It is important to consult with your healthcare professional before using this device, especially if you have any medical conditions or concerns. This device is intended to supplement, not replace, professional medical advice.

Disclaimer:

The manufacturer is not responsible for any misuse or improper use of the device. Use of the device without consulting a doctor may result in inaccurate readings or other adverse effects. Always rely on professional medical advice and do not rely solely on the device.

Contact Information:

For technical support or assistance, please contact our customer service team at **ketosmartapp@gmail.com**.